


# POSTURAL ACCESSORIES



- Anti-slip & waterproof seat cushion made of multi-density foams resting on a firm base for an optimal dispersion of the pelvic pressure zones
- A unique combination of fabrics creates a true anti-slip effect
- Makes it more difficult for the user to slide off of a wheelchair.
- It helps retain the proper position of the pelvis and to maintain the desired alignment of the spine
- Outer & inner layers of waterproof barrier for added protection
- Available in Child & Adult configurations with 18 dimensions
- Compatible with postural backrests : FLEXUS™, BACKUS™ & CURVUS™
- Hand washable cover


SPECIFICATIONS			DEPTH*					
16 x 18"		WIDTH**	11 to 12"	13 to 14"	15 to 16"	17 to 18"	19 to 20"	
Weight of Cushion	2 lb / 1.4kg	11 to 12"	STOP111212	STOP111214	STOP111216			
		13 to 14"	STOP131412	STOP131414	STOP131416			
Thickness	3 3/8" / 9cm	15 to 16"			STOP151616	STOP151618	STOP151620	
		17 to 18"			STOP171816	STOP171818	STOP171820	
Weight limitation	200 lbs / 120kg	19 to 20"			STOP192016	STOP192018	STOP192020	
		21 to 22"			STOP212216	STOP212218	STOP212220	

\* The actual depths are 12-14-16-18-20", but the rounded edges allow to accommodate a shorter depth of 1".



\*\* The actual widths are 11 1/2 - 13 1/2 - 15 1/2 - 17 1/2 - 19 1/2 - 21 1/2", and each cushion is intended for two widths of wheelchair.

N.B. (1) THE INSTALLATION IS A FULL WIDTH FLAT BASE



(2) TO ORDER AN OUTER CUSHION COVER, PLEASE ADD THE LETTER < C > TO THE PRODUCT CODE.

(3) THE VELCRO™ OF THE COVER HAS A LOOP FINISH.

**Posi**  
**tech**  
**innovation™**

100-1566 Rue Nationale,  
Terrebonne, Québec,  
Canada, J6W 0E2

TEL: 1 450-964-3232

1 877-964-3232

www.positechinnovation.com  
service@positechinnovation.com