POSTURAL ACCESSORIES







BREASTPLATES WITH ANTI-ELASTIC

- Backrest fastening system to secure the proper posture of the chest while sitting in a wheelchair
- The anti-elastic feature keeps the wheelchair user in a more secure position
- Adjustable in length and width
- Available in 5 models:3 sizes & 2 shapes:'H' and 'X'
- Modular system using compatible postural backrests: FLEXUS™, BACKUS™ & CURVUS™
- Very soft and hand washable

SPECIFICATIONS		'H'	′X′	
LENGTH	Small - 12" (30cm)	STAR1182		
	Medium - 15" (38cm)	STAR1183		
	Large - 18" (46cm)	STAR1184		
	Medium - 13" (33cm)		STAR1193	
	Large - 16" (41cm)		STAR1194	
		HARDWARE		
Breastplate anchoring system with mobile cursors		STAR1151		
Backrest fastener		SUMM1110		
Reinforcement bar		SUMM1106		
Mobile cursors to anchor breastplate		STAR1152		





POSTURE BELTS

- A versatile 4-anchor-point-buckle belt that procures, depending how it is positioned onto the wheelchair, additional support for an optimal posture of the back or the pelvis
- Available in 3 sizes
- Modular system using cooperable postural back estar FLD/USTM, BACKUSTM & CURVUSTM
- Very soft and hand washable

SPECIFICATIONS		
PELVIC WIDTH	Small - 10 to 13" (25 - 33cm)	STAR1162
	Medium - 14 to 17" (35 - 43cm)	STAR1163
	Large - 18 to 21" (45 - 53cm)	STAR1164



CORPORATE HEADQUARTERS

2-1096 rue Lévis, Terrebonne, Québec, Canada J6W4L1

TEL: 1 450-964-3232 1 877-964-3232 FAX: 1 450-964-3378

www.positechinnovation.com info@positechinnovation.com





. Produc Caaog 20

- In order to optimize the comfort and support, FORMUS™ is a back-support foam lumbar insert
- With FLEXUS[™], FORMUS[™] is inserted in between the rigid structure and the upholstered cushion, and secure on the Velcro strip
- With BACKUSTM & CURVUSTM, FORMUSTM is inserted inside the upholstered cushion on the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the patient's back side and secure to the declaration of the declaration of the patient's back side and secure to the declaration of the declaration of
- Available in 2 profiled shapes: contoured and bevele@DULT
- Each profiled shape is available in 4 sizes
- Modular system using compatible postural backrests: FLEXUS™, BACKUS™ & CURVUS™
- Wipe with damp cloth to clean

SPECIFICATIONS	CONTOURED		BEVELED	
SI ECITICATIONS	ADULT	CHILD	ADULT	CHILD
1" - Small - 6 x 12" (20 x 30cm)	FORM1240			
1 1/2" - Medium - 6 x 12" (20 x 30cm)	FORM1241			
1/2" - X-small - 5 x 8 " (12 x 20cm)		FORM1270		
1" - Small - 5 x 8" (12 x 20 cm)		FORM1271		
1" - Small - 6 x 12" (20 x 30cm)			CAR41863 2	
1 1/2" - Medium - 6 x 12" (20 x 30cm)			FORM1243	
1/2" - X-small - 5 x 8" (12 x 20 cm)				FORM1272
1" - Small - 5 x 8" (12 x 20 cm)				FORM1273